



## Gloves

Gloves are one of the top requested items on everyone's safety lists, we have options and sizes to fit the smallest to largest hands.

- Gloves should fit snug, but comfortable on the hand
- Gloves should be peeled off from the wrist down and off the fingers, so that they are inside out by the time they are off the hand

Flinn has a full selection of gloves available to suit every need from providing simple hand coverage to using while cleaning and disinfecting classrooms and more. As with any personal protective equipment, it is important to select gloves that are appropriate for the chemical hazards. One type of glove does NOT protect against all types of chemicals!

Chemical resistance of gloves is measured in three ways: degradation, breakthrough time, and permeation rate. For heavy-duty or longterm use you must consider not only the glove material but also the thickness, size and fit of the gloves. Review our [Glove Selection Guide](#) for details!

### Recommended Products:

	PRODUCT NAME
AP7080*	Disposable Powder-Free Nitrile Gloves Pkg. of 100. Sizes XL-XS
AP3216*	Disposable Polyethylene Gloves Pkg. of 100. Sizes XI-S
SE1032	Disposable Polyethylene Gloves Pkg. of 500
AP3232	Playtex-Type Gloves 1 Pair
AP4429*	Latex Gloves Pkg. of 100. Sizes XI-S
AP4828*	Powder-Free Latex Glove Pkg. of 100. Sizes XI-S

For a Customized Solution for your school or district, contact a Flinn sales representatives today at 1-800-452-1261

\*each size has it's own product number; numbers listed are for large sizes

### WHEN TO WEAR GLOVES:

The [CDC recommends](#) using gloves when cleaning and when caring for someone who may potentially be sick.

Schools should provide gloves to any staff or students who may come in contact with a person showing signs of illness, or when taking daily temperature checks.

And as always, students participating in Lab Experiments should wear gloves during all science activities.

